

Ultimate Comfort SMART

DENTAL TRADITION
130
YEARS

Ritter®
[THE DENTAL EXPERTS]

new



H-Version

hanging version

- up to 5 instruments
- touchpanel with x-ray viewer
- 3 programmable patient chair positions



S-Version

swinging version

- up to 5 instruments
- touchpanel with x-ray viewer
- 3 programmable patient chair positions

Assistant Element

- Integrated touch panel
- Up to 4 instruments
- Vacuum package small saliva ejector/high volume ejector
- 3F syringe
- Moveable arm



Please note: Pictures show equipment which might be available only as an option.
Differences between samples and original surface are technically imposed.
Alterations due to new technical developments reserved. V17.00.00

Ultimate Comfort SMART

new

DENTAL TRADITION
130
YEARS

Ritter®
[THE DENTAL EXPERTS]

Features and options

● STANDARD ○ OPTION — NOT AVAILABLE

Patient chair

Trendelenburg movement
Storage of positions
Safety switch
2 articulated headrest mechanically
Left and right armrest
Seamless upholstery
Multifunctional control

Dentist element

3 Functional syringe
3 outlets for turbine or air motor without light
Fibre optic outlets
Motor Infinity brushless LED
ZEG Kolibri S or E
Small x-ray film viewer
Recording for 2 standard trays

Assistant element

Operating panel
Moveable assistant element
3 Functional syringe
Large suction hose
Small suction hose
LED polymerization handpiece

Water unit

Moveable ceramic cuspidor
ABS Aqua Bottle System
Heating for cup filler
Connection for external suction
Place selection valve
Venturi System

Operating light

Ritter PlanetLite Hologen
Ritter CoolLite LED

Stools

Mobilo Flex 50
Mobilo Flex 50 with footring

Multi media

Intraoral camera

Foot controller

Modern design for chair movement and instruments:

- Large foot bow for easy repositioning
- Manual control of chair movements
- Up/down/back and forth

Program control of chair movements:

- Spitting position/last treatment position
- Free programmable chair positions
- Exit position

Control of instruments:

- Speed/power control for the active instrument

Buttons:

- Spray water at instrument
- Chip blower function
- Optional foot controller pedal type



Stool

- Backrest with fixed position or dynamic movement
- Ergonomic and comfortable seating
- Symmetrical backrest
- Gas lift cylinder
- Chrome legs
- Five dual castors for stability
- Colour-matching upholstery available in multiple colours
- Available with or without footring

Upholstery colours:

Multiple fresh and modern colors are available.

